




Dablood (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)

Clinical Trials Evaluation

This clinical trial evaluation unit is the municipal public health center of Jiangsu Province and the affiliated hospital of Sun Yat-sen University. The entire clinical trial process is carried out under strict control and are conducted by professional trained testers





Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit

Clinical Trials Evaluation

Sampling Distribution

This trial used a blinded and controlled clinical trial design to conduct a clinical trial of the "COVID-19 (SARS-CoV-2) IgG / IgM antibody detection kit (colloidal gold method)" produced by Jiangsu Dablood Pharmaceutical Co., Ltd. The clinical trial unit is a public health center of a city in Jiangsu Province and an affiliated hospital of Sun Yat-sen University. The entire clinical trial process is carried out under strict control, and detection and analysis are conducted by professional trained testers. This clinical trial included 424 cases, including 424 serum specimens. Among the 424 cases, 183 were male (43.16%), 232 were female (54.72%), and 9 were missing in gender information (2.12%); 43 patients (10.14%) with age $0 < Y \leq 20$, 20 < $Y \leq 40$ 112 cases (28.06%) for ≤ 40 , 193 cases (47.64%) for $40 < Y \leq 60$, 51 cases (12.03%) for $60 < Y \leq 80$, 9 cases (2.12%) for $Y > 80$. There were 16 unknown age cases (8.92%); 63 of the 424 cases were diagnosed as negative cases, and 311 cases were tested positive and clinically diagnosed positive. 36 cases physical examination in good health, 2 respiratory syncytial virus (RSV) cases, 5 influenza A viruses, 5 (Flu A) cases, 5 influenza B virus (Flu B) cases, 1 cancer tumor, and other diseases 1 case. The test results show that the positive case-control group and the test group have a good consistency, and at the same time, the negative case did not detect a positive result with both reagents.



Appendix 1.

**Jiangsu Provincial Public Health Center
positive sample clinical data table**

Sample No.	Clinical diagnosis	Control reagent IgG results	Control reagent IgM results	Dablood reagent IgG results	Dablood reagent IgM results
1	SARS-CoV-2	+	-	+	-
2	SARS-CoV-2	+	-	+	-
3	SARS-CoV-2	+	-	+	-
4	SARS-CoV-2	+	-	-	-
5	SARS-CoV-2	+	-	+	-
6	SARS-CoV-2	+	-	+	-
7	SARS-CoV-2	+	-	-	-
8	SARS-CoV-2	+	-	+	-
9	SARS-CoV-2	+	-	+	-
10	SARS-CoV-2	+	-	-	-
11	SARS-CoV-2	+	-	+	-
12	SARS-CoV-2	+	-	+	-
13	SARS-CoV-2	+	-	+	-
14	SARS-CoV-2	+	-	+	-
15	SARS-CoV-2	+	-	+	-
16	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



17	SARS-CoV-2	+	-	-	-
18	SARS-CoV-2	+	-	+	-
19	SARS-CoV-2	+	-	+	-
20	SARS-CoV-2	+	-	+	-
21	SARS-CoV-2	+	-	+	-
22	SARS-CoV-2	+	-	+	-
23	SARS-CoV-2	+	-	+	-
24	SARS-CoV-2	+	-	+	-
25	SARS-CoV-2	+	-	+	-
26	SARS-CoV-2	+	-	+	-
27	SARS-CoV-2	+	-	+	-
28	SARS-CoV-2	+	-	+	-
29	SARS-CoV-2	+	-	+	-
30	SARS-CoV-2	+	-	+	-
31	SARS-CoV-2	+	-	+	-
32	SARS-CoV-2	+	-	-	-
33	SARS-CoV-2	+	-	+	-
34	SARS-CoV-2	+	-	+	-
35	SARS-CoV-2	+	-	+	-
36	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



37	SARS-CoV-2	+	-	+	-
38	SARS-CoV-2	+	-	+	-
39	SARS-CoV-2	+	-	-	-
40	SARS-CoV-2	+	-	+	-
41	SARS-CoV-2	+	-	+	-
42	SARS-CoV-2	+	-	+	-
43	SARS-CoV-2	+	-	+	-
44	SARS-CoV-2	+	-	+	-
45	SARS-CoV-2	+	-	+	-
46	SARS-CoV-2	+	-	+	-
47	SARS-CoV-2	+	-	+	-
48	SARS-CoV-2	+	-	+	-
49	SARS-CoV-2	+	-	+	-
50	SARS-CoV-2	+	-	+	-
51	SARS-CoV-2	+	-	+	-
52	SARS-CoV-2	+	-	+	-
53	SARS-CoV-2	+	-	+	-
54	SARS-CoV-2	+	-	+	-
55	SARS-CoV-2	+	-	+	-
56	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



57	SARS-CoV-2	+	-	+	-
58	SARS-CoV-2	+	-	+	-
59	SARS-CoV-2	+	-	+	-
60	SARS-CoV-2	+	-	+	-
61	SARS-CoV-2	+	-	-	-
62	SARS-CoV-2	+	-	+	-
63	SARS-CoV-2	+	-	+	-
64	SARS-CoV-2	+	-	+	-
65	SARS-CoV-2	+	-	+	-
66	SARS-CoV-2	+	-	+	-
67	SARS-CoV-2	+	-	+	-
68	SARS-CoV-2	+	-	+	-
69	SARS-CoV-2	+	-	+	-
70	SARS-CoV-2	+	-	+	-
71	SARS-CoV-2	+	-	+	-
72	SARS-CoV-2	+	-	+	-
73	SARS-CoV-2	+	-	+	-
74	SARS-CoV-2	+	-	-	-
75	SARS-CoV-2	+	-	+	-
76	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



77	SARS-CoV-2	+	-	+	-
78	SARS-CoV-2	+	-	+	-
79	SARS-CoV-2	+	-	+	-
80	SARS-CoV-2	+	-	+	-
81	SARS-CoV-2	+	-	+	-
82	SARS-CoV-2	+	-	+	-
83	SARS-CoV-2	+	-	+	-
84	SARS-CoV-2	+	-	+	-
85	SARS-CoV-2	+	-	+	-
86	SARS-CoV-2	+	-	+	-
87	SARS-CoV-2	+	-	+	-
88	SARS-CoV-2	+	-	+	-
89	SARS-CoV-2	+	-	+	-
90	SARS-CoV-2	+	-	+	-
91	SARS-CoV-2	+	-	+	-
92	SARS-CoV-2	+	-	+	-
93	SARS-CoV-2	+	-	-	-
94	SARS-CoV-2	+	-	+	-
95	SARS-CoV-2	+	-	+	-
96	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



97	SARS-CoV-2	+	-	+	-
98	SARS-CoV-2	+	-	+	-
99	SARS-CoV-2	+	-	+	-
100	SARS-CoV-2	+	-	+	-
101	SARS-CoV-2	+	-	+	-
102	SARS-CoV-2	+	-	+	-
103	SARS-CoV-2	+	-	+	-
104	SARS-CoV-2	+	-	+	-
105	SARS-CoV-2	+	-	+	-
106	SARS-CoV-2	+	-	+	-
107	SARS-CoV-2	+	-	+	-
108	SARS-CoV-2	+	-	+	+
109	SARS-CoV-2	+	-	+	-
110	SARS-CoV-2	+	-	+	-
111	SARS-CoV-2	+	-	+	-
112	SARS-CoV-2	+	-	+	-
113	SARS-CoV-2	+	-	+	-
114	SARS-CoV-2	+	-	+	-
115	SARS-CoV-2	+	-	+	-
116	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



117	SARS-CoV-2	+	-	+	-
118	SARS-CoV-2	+	-	+	-
119	SARS-CoV-2	-	+	-	+
120	SARS-CoV-2	-	+	-	+
121	SARS-CoV-2	-	+	-	+
122	SARS-CoV-2	-	+	+	+
123	SARS-CoV-2	-	+	-	+
124	SARS-CoV-2	-	+	-	+
125	SARS-CoV-2	-	+	-	+
126	SARS-CoV-2	-	+	+	+
127	SARS-CoV-2	-	+	-	-
128	SARS-CoV-2	-	+	-	+
129	SARS-CoV-2	-	+	-	+
130	SARS-CoV-2	-	+	-	+
131	SARS-CoV-2	-	+	-	+
132	SARS-CoV-2	-	+	-	+
133	SARS-CoV-2	-	+	-	+
134	SARS-CoV-2	-	+	-	+
135	SARS-CoV-2	-	+	-	-
136	SARS-CoV-2	-	+	-	+

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



137	SARS-CoV-2	-	+	-	+
138	SARS-CoV-2	-	+	-	+
139	SARS-CoV-2	-	+	-	+
140	SARS-CoV-2	+	+	-	+
141	SARS-CoV-2	+	+	+	+
142	SARS-CoV-2	+	+	+	+
143	SARS-CoV-2	+	+	-	+
144	SARS-CoV-2	+	+	+	+
145	SARS-CoV-2	+	+	+	+
146	SARS-CoV-2	+	+	+	+
147	SARS-CoV-2	+	+	+	-
148	SARS-CoV-2	+	+	+	+
149	SARS-CoV-2	+	+	+	+
150	SARS-CoV-2	+	+	+	+
151	SARS-CoV-2	-	-	-	-
152	SARS-CoV-2	-	-	-	-
153	SARS-CoV-2	-	-	-	-
154	SARS-CoV-2	-	-	-	-
155	SARS-CoV-2	-	-	-	-
156	SARS-CoV-2	-	-	-	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



157	SARS-CoV-2	-	-	-	-
158	SARS-CoV-2	-	-	+	-
159	SARS-CoV-2	-	-	-	-
160	SARS-CoV-2	-	-	-	-
161	SARS-CoV-2	-	-	-	-
162	SARS-CoV-2	-	-	-	-
163	SARS-CoV-2	-	-	-	-
164	SARS-CoV-2	-	-	-	-
165	SARS-CoV-2	-	-	-	-
166	SARS-CoV-2	-	-	-	-
167	SARS-CoV-2	-	-	-	-
Total		Positive : 129 Negative: 18	Positive: 32 Negative: 135	Positive: 121 Negative: 46	Positive: 30 Negative: 137



Appendix 2.

Sun Yat-Sen University Affiliated Hospital Positive Samples Clinical Data Table

Sample No.	Clinical diagnosis	Control reagent IgG results	Control reagent IgM results	Dablood reagent IgG results	Dablood reagent IgM results
1	SARS-CoV-2	+	-	+	-
2	SARS-CoV-2	+	-	+	-
3	SARS-CoV-2	+	-	+	-
4	SARS-CoV-2	+	-	+	-
5	SARS-CoV-2	+	-	+	-
6	SARS-CoV-2	+	-	+	-
7	SARS-CoV-2	+	-	+	-
8	SARS-CoV-2	+	-	-	-
9	SARS-CoV-2	+	-	+	-
10	SARS-CoV-2	+	-	+	-
11	SARS-CoV-2	+	-	-	-
12	SARS-CoV-2	+	-	+	-
13	SARS-CoV-2	+	-	+	-
14	SARS-CoV-2	+	-	+	-
15	SARS-CoV-2	+	-	+	-
16	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



17	SARS-CoV-2	+	-	+	-
18	SARS-CoV-2	+	-	-	-
19	SARS-CoV-2	+	-	+	-
20	SARS-CoV-2	+	-	+	-
21	SARS-CoV-2	+	-	+	-
22	SARS-CoV-2	+	-	+	-
23	SARS-CoV-2	+	-	+	-
24	SARS-CoV-2	+	-	+	-
25	SARS-CoV-2	+	-	+	-
26	SARS-CoV-2	+	-	-	-
27	SARS-CoV-2	+	-	+	-
28	SARS-CoV-2	+	-	+	-
29	SARS-CoV-2	+	-	+	-
30	SARS-CoV-2	+	-	+	-
31	SARS-CoV-2	+	-	+	-
32	SARS-CoV-2	+	-	+	-
33	SARS-CoV-2	+	-	+	-
34	SARS-CoV-2	+	-	+	-
35	SARS-CoV-2	+	-	-	-
36	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



37	SARS-CoV-2	+	-	+	-
38	SARS-CoV-2	+	-	+	-
39	SARS-CoV-2	+	-	-	-
40	SARS-CoV-2	+	-	+	-
41	SARS-CoV-2	+	-	+	-
42	SARS-CoV-2	+	-	+	-
43	SARS-CoV-2	+	-	+	-
44	SARS-CoV-2	+	-	+	-
45	SARS-CoV-2	+	-	+	-
46	SARS-CoV-2	+	-	+	-
47	SARS-CoV-2	+	-	+	-
48	SARS-CoV-2	+	-	+	-
49	SARS-CoV-2	+	-	+	-
50	SARS-CoV-2	+	-	+	-
51	SARS-CoV-2	+	-	+	-
52	SARS-CoV-2	+	-	+	-
53	SARS-CoV-2	+	-	+	-
54	SARS-CoV-2	+	-	+	-
55	SARS-CoV-2	+	-	+	-
56	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



57	SARS-CoV-2	+	-	+	-
58	SARS-CoV-2	+	-	-	-
59	SARS-CoV-2	+	-	+	-
60	SARS-CoV-2	+	-	+	-
61	SARS-CoV-2	+	-	+	-
62	SARS-CoV-2	+	-	+	-
63	SARS-CoV-2	+	-	+	-
64	SARS-CoV-2	+	-	+	-
65	SARS-CoV-2	+	-	+	-
66	SARS-CoV-2	+	-	+	-
67	SARS-CoV-2	+	-	+	-
68	SARS-CoV-2	+	-	+	-
69	SARS-CoV-2	+	-	+	-
70	SARS-CoV-2	+	-	+	-
71	SARS-CoV-2	+	-	+	-
72	SARS-CoV-2	+	-	+	-
73	SARS-CoV-2	+	-	+	-
74	SARS-CoV-2	+	-	+	-
75	SARS-CoV-2	+	-	+	-
76	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



77	SARS-CoV-2	+	-	+	-
78	SARS-CoV-2	+	-	+	-
79	SARS-CoV-2	+	-	+	-
80	SARS-CoV-2	+	-	+	-
81	SARS-CoV-2	+	-	+	-
82	SARS-CoV-2	+	-	-	+
83	SARS-CoV-2	+	-	+	-
84	SARS-CoV-2	+	-	+	-
85	SARS-CoV-2	+	-	+	-
86	SARS-CoV-2	+	-	+	-
87	SARS-CoV-2	+	-	+	-
88	SARS-CoV-2	+	-	+	-
89	SARS-CoV-2	+	-	-	-
90	SARS-CoV-2	+	-	+	-
91	SARS-CoV-2	+	-	+	-
92	SARS-CoV-2	+	-	+	-
93	SARS-CoV-2	+	-	+	-
94	SARS-CoV-2	+	-	+	-
95	SARS-CoV-2	+	-	+	-
96	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



97	SARS-CoV-2	+	-	+	-
98	SARS-CoV-2	+	-	+	-
99	SARS-CoV-2	+	-	+	-
100	SARS-CoV-2	+	-	+	+
101	SARS-CoV-2	+	-	+	-
102	SARS-CoV-2	+	-	+	-
103	SARS-CoV-2	+	-	+	-
104	SARS-CoV-2	+	-	+	-
105	SARS-CoV-2	+	-	+	-
106	SARS-CoV-2	+	-	+	-
107	SARS-CoV-2	+	-	+	-
108	SARS-CoV-2	+	-	+	-
109	SARS-CoV-2	+	-	+	-
110	SARS-CoV-2	+	-	+	-
111	SARS-CoV-2	+	-	+	-
112	SARS-CoV-2	+	-	+	-
113	SARS-CoV-2	+	-	+	-
114	SARS-CoV-2	+	-	+	-
115	SARS-CoV-2	+	-	+	-
116	SARS-CoV-2	+	-	+	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



117	SARS-CoV-2	+	-	+	-
118	SARS-CoV-2	+	-	+	-
119	SARS-CoV-2	-	+	-	+
120	SARS-CoV-2	-	+	-	+
121	SARS-CoV-2	-	+	-	+
122	SARS-CoV-2	-	+	+	+
123	SARS-CoV-2	-	+	-	+
124	SARS-CoV-2	-	+	-	+
125	SARS-CoV-2	-	+	-	-
126	SARS-CoV-2	-	+	-	+
127	SARS-CoV-2	-	+	-	+
128	SARS-CoV-2	-	+	+	+
129	SARS-CoV-2	-	+	-	+
130	SARS-CoV-2	-	+	-	+
131	SARS-CoV-2	-	+	-	+
132	SARS-CoV-2	+	+	-	+
133	SARS-CoV-2	+	+	+	+
134	SARS-CoV-2	+	+	+	+
135	SARS-CoV-2	+	+	+	-
136	SARS-CoV-2	-	-	-	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



137	SARS-CoV-2	-	-	+	-
138	SARS-CoV-2	-	-	-	-
139	SARS-CoV-2	-	-	-	-
140	SARS-CoV-2	-	-	-	-
141	SARS-CoV-2	-	-	-	-
142	SARS-CoV-2	-	-	-	-
143	SARS-CoV-2	-	-	-	-
144	SARS-CoV-2	-	-	-	-
Total		Positive: 122 Negative: 22	Positive: 17 Negative: 127	Positive: 115 Negative: 29	Positive: 17 Negative: 127



Appendix 3.

Nanjing Gulou Hospital Negative Samples Clinical Data Table

Sample No.	Clinical diagnosis	Control reagent IgG results	Control reagent IgM results	Dablood reagent IgG results	Dablood reagent IgM results
1	Mycoplasma pneumonia	-	-	-	-
2	Common cold, fever	-	-	-	-
3	Streptococcus pneumonia	-	-	-	-
4	Chlamydia pneumonia	-	-	-	-
5	Bronchiolitis	-	-	-	-
6	Bronchopneumonia	-	-	-	-
7	Common cold, fever	-	-	-	-
8	Chlamydia pneumonia	-	-	-	-
9	Bronchopneumonia	-	-	-	-
10	Common cold, fever	-	-	-	-
11	Chlamydia pneumonia	-	-	-	-
12	Common cold, fever	-	-	-	-



13	Mycoplasma pneumonia	-	-	-	-
14	Chlamydia pneumonia	-	-	-	-
15	Bronchiolitis	-	-	-	-
16	Chlamydia pneumonia	-	-	-	-
17	Bronchiolitis	-	-	-	-
18	Common cold, fever	-	-	-	-
19	Interstitial pneumonia	-	-	-	-
20	Common cold, fever	-	-	-	-
21	Interstitial pneumonia	-	-	-	-
22	Bronchiolitis	-	-	-	-
23	Mycoplasma pneumonia	-	-	-	-
24	Chlamydia pneumonia	-	-	-	-
25	Bronchopneumonia	-	-	-	-
26	Common cold, fever	-	-	-	-
27	Bronchiolitis	-	-	-	-
28	Common cold, fever	-	-	-	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



29	Bronchiolitis	-	-	-	-
30	Mycoplasma pneumonia	-	-	-	-
31	Common cold, fever	-	-	-	-
32	Chlamydia pneumonia	-	-	-	-
33	Common cold, fever	-	-	-	-
34	Bronchiolitis	-	-	-	-
35	Bronchiolitis	-	-	-	-
36	Interstitial pneumonia	-	-	-	-
37	Bronchiolitis	-	-	-	-
38	Chlamydia pneumonia	-	-	-	-
39	Common cold, fever	-	-	-	-
40	Bronchopneumonia	-	-	-	-
41	Mycoplasma pneumonia	-	-	-	-
42	Common cold, fever	-	-	-	-
43	Bronchiolitis	-	-	-	-
44	Mycoplasma pneumonia	-	-	-	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



45	Bronchiolitis	-	-	-	-
46	Chlamydia pneumonia	-	-	-	-
47	Interstitial pneumonia	-	-	-	-
48	Bronchiolitis	-	-	-	-
49	Common cold, fever	-	-	-	-
50	Bronchopneumonia	-	-	-	-
51	Common cold, fever	-	-	-	-
52	Mycoplasma pneumonia	-	-	-	-
53	Bronchopneumonia	-	-	-	-
54	Common cold, fever	-	-	-	-
55	Bronchopneumonia	-	-	-	-
56	Mycoplasma pneumonia	-	-	-	-
57	Interstitial pneumonia	-	-	-	-
58	Common cold, fever	-	-	-	-
59	Chlamydia pneumonia	-	-	-	-



60	Interstitial pneumonia	-	-	-	-
61	Common cold, fever	-	-	-	-
62	Chlamydia pneumonia	-	-	-	-
63	Common cold, fever	-	-	-	-
64	Healthy case	-	-	-	-
65	Healthy case	-	-	-	-
66	Healthy case	-	-	-	-
67	Healthy case	-	-	-	-
68	Healthy case	-	-	-	-
69	Healthy case	-	-	-	-
70	Healthy case	-	-	-	-
71	Healthy case	-	-	-	-
72	Healthy case	-	-	-	-
73	Healthy case	-	-	-	-
74	Healthy case	-	-	-	-
75	Healthy case	-	-	-	-
76	Healthy case	-	-	-	-
77	Healthy case	-	-	-	-
78	Healthy case	-	-	-	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



79	Healthy case	-	-	-	-
80	Healthy case	-	-	-	-
81	Healthy case	-	-	-	-
82	Healthy case	-	-	-	-
83	Healthy case	-	-	-	-
84	Healthy case	-	-	-	-
85	Healthy case	-	-	-	-
86	Healthy case	-	-	-	-
87	Healthy case	-	-	-	-
88	Healthy case	-	-	-	-
89	Healthy case	-	-	-	-
90	Healthy case	-	-	-	-
91	Healthy case	-	-	-	-
92	Healthy case	-	-	-	-
93	Healthy case	-	-	-	-
94	Healthy case	-	-	-	-
95	Healthy case	-	-	-	-
96	Healthy case	-	-	-	-
97	Healthy case	-	-	-	-
98	Healthy case	-	-	-	-

Dablood Pharmaceutical (SARS-CoV-2) IgG / IgM test kit (Colloidal Gold)



99	Healthy case	-	-	-	-
100	Respiratory syncytial virus	-	-	-	-
101	Respiratory syncytial virus	-	-	-	-
102	Influenza A virus	-	-	-	-
103	Influenza A virus	-	-	-	-
104	Influenza A virus	-	-	-	-
105	Influenza A virus	-	-	-	-
106	Influenza A virus	-	-	-	-
107	Influenza B virus	-	-	-	-
108	Influenza B virus	-	-	-	-
109	Influenza B virus	-	-	-	-
110	Influenza B virus	-	-	-	-
111	Influenza B virus	-	-	-	-
112	Lung cancer case	-	-	-	-
113	TYPE I DIABETES	-	-	-	-
Total		Positive: 0 Negative: 113	Positive: 0 Negative: 113	Positive: 0 Negative: 113	Positive: 0 Negative: 113



3 IgG Statistical description

testing facility		IgG Test		Total			
		+	-				
Jiangsu Province	Dablood Test	+	Subtotal	118	3	121	
		-	Subtotal	11	86	97	
Public Health Center		Total		Subtotal	129	89	218
Sun Yat-Sen University Affiliated Hospital	Dablood Test	+	Subtotal	113	3	116	
		-	Subtotal	9	81	90	
		Total		Subtotal	122	84	206
Total	Dablood Test	+	Subtotal	232	6	238	
		-	Subtotal	19	167	186	
		Total		Subtotal	251	173	424

4 Coincidence Rate

+ Coincidence Rate : $232 / 251 \times 100\% = 92.43\%$

- Coincidence Rate : $167 / 173 \times 100\% = 96.53\%$

Overall Coincidence Rate : $(232 + 167) / 424 \times 100\% = 94.10\%$



5 IgM Statistical description

testing facility				IgM Test		Total	
		+	-	+	-		
Jiangsu Province	Dablood Test	+	Subtotal	29	1	30	
		-	Subtotal	3	185	188	
Public Health Center		Total		Subtotal	32	186	218
Sun Yat-Sen University Affiliated Hospital	Dablood Test	+	Subtotal	15	2	17	
		-	Subtotal	2	187	189	
Total		Total		Subtotal	17	189	206
Total	Dablood Test	+	Subtotal	44	3	47	
		-	Subtotal	5	372	377	
Total		Total		Subtotal	49	375	424

6 Coincidence Rate

+ Coincidence Rate : $44 / 47 \times 100\% = 93.62\%$

- Coincidence Rate : $372 / 377 \times 100\% = 98.67\%$

Overall Coincidence Rate : $(44 + 372) / 424 \times 100\% = 98.11\%$